

# Dosing Package pH Meter Package Conductivity Package



Complete packages for  
professional liquid handling and  
pH and conductivity measurement

## Dosing Package for liquid handling ...

02

The Dosing Package is the practical and proven solution for precise, reliable liquid handling in accordance with directives, both in regulated and non-regulated environments.

You can rely on the Dosing Package for:

- Standard applications such as dosing, diluting, taking aliquots, pipetting
- The preparation of standard solutions from any given sample size
- The preparation of mixtures from various stock solutions



### The Dosing Package includes

- 1 x 846 Dosing Interface
- 1 x 800 Dosing Unit
- 1 x 900 Touch Control
- 1 x 20 mL dosing unit

### Your benefit

- Fulfills requirements of FDA regulation 21 CFR Part 11
- Each application can be programmed as a method and started at the press of a button
- Complete monitoring and traceability of applications and results
- Thanks to USB and Ethernet interface, a wide range of options for data backup (LIMS, PDF report, USB flash drive, local or network printer)
- Completely compatible with other Metrohm Titrandos instruments (pH and conductivity measurement, titration, automation)

### Order number

- Dosing Package: 2.846.1010

# Dosing Package


## ... and pH Meter Package for pH-value/ ion measurement and Conductivity Package for conductivity measurement

Both the 867 pH Module and the 856 Conductivity Module with Touch Control are measuring instruments that meet the highest standards.

03

They are particularly suitable for:

- Measurement of the pH value, the ORP, the temperature or an ion concentration
- pH measurement in accordance with USP 791
- Ion measurement using fully automated standard addition (with optional 800 Dosino)
- Conductivity measurement in accordance with USP <645>



### Your benefit

- Fulfills requirements pursuant to FDA regulation 21 CFR Part 11
- Each application can be programmed as a method and started at the press of a button
- Complete monitoring and traceability of applications and results
- Thanks to USB and Ethernet interface, a wide range of opportunities for data backup (LIMS, PDF report, USB flash drive, local or network printer)
- Completely compatible with other Metrohm Titrandos instruments (titration, automation)

An optional 800 Dosino provides the 867 pH Module and the 856 Conductivity Module the same liquid handling functionality as the 846 Dosing Interface

### Order numbers

- **pH Meter Package: 2.867.1410**
- **Conductivity Package: 2.856.0120**

# Ordering options

## 2.846.1010 Dosing Package

- 846 Dosing Interface
- 800 Dosino
- 900 Touch Control
- 6.2061.010 Bottle holder for Dosino
- 6.2065.000 Stacking frame for 846
- 6.3032.220 Dosing unit 20 mL
- 6.6064.010 USB flash drive for 900

## 2.867.1410 pH Meter Package

- 867 pH Module
- 900 Touch Control
- 6.6064.010 USB flash drive for 900  
including stand and electrode holder (sensor not included\*)

## 2.856.0120 Conductivity Module with Touch Control and steel cell

- Conductivity Module
- 900 Touch Control
- 6.0916.040 Conductivity measuring cell, stainless steel with  
Pt1000 m C=0.1
- 6.6064.010 USB flash drive for 900  
including stand and electrode holder

---

### Optional accessories

- 2.800.0010 800 Dosino
- 2.801.0040 Magnetic stirrer MSB
- 6.3032.120 Dosing unit 2 mL
- 6.3032.150 Dosing unit 5 mL
- 6.3032.210 Dosing unit 10 mL
- 6.3032.220 Dosing unit 20 mL
- 6.3032.250 Dosing unit 50 mL

---

\*You'll find the right electrode for your application under  
[www.metrohm.com/com/products/electrodes](http://www.metrohm.com/com/products/electrodes)  
or with our Electrode Finder App!



[www.metrohm.com](http://www.metrohm.com)

